

Tech Laser Therapy 302



LASER THERAPY UNIT

Anti inflammatory



Wound healing



Pain management

Tech Laser Therapy 302

LASER THERAPY UNIT

TECH LASER THERAPY 302 is a advanced Microprocessor based Laser Therapy Equipment. It will provide top therapeutic efficacy, for your satisfaction and for your patient's comfort.

TECH LASER THERAPY 302 is a Portable user-friendly unit with zoomed LCD panel with backlights that displays the TREATMENT TIME and FREQUENCY. It also has an attractive BARGRAPH display and BUZZER INDICATION for acupuncture detection. TECH PORTABLE LASER 302 provides two types of PROBES, VISIBLE RED LASER AND INFRARED LASER.

FEATURES

- ◆ Computerised LCD screen with controls gives messages / information to the user on various adjustments.
- ◆ Two Probes - Visible Red Probe and Infra Red Probe
- ◆ Automatic electronic treatment timer.
- ◆ Pulsed and continuous output.
- ◆ Frequency 0.6 to 300 Hz
- ◆ Moulded ABS Body Portable Table Top

TECHNICAL DATA

Mains supply	:	230V/AC/50 Hz
Power Consumption	:	15 Watts
Mains fuse	:	1 Amp
Treat Time	:	max-99 min
Frequency	:	0.6 to 300 Hz
Wavelength Infrared	:	980 nm / 100 mw
Visible Red Laser	:	655 nm / 20 mw