

# TECH LASER THERAPY SS



## LASER THERAPY UNIT

Anti inflammatory



Wound healing



Pain management



# TECH LASER THERAPY SS

## LASER THERAPY UNIT

TECH Laser is a versatile, Solid State Laser Output with visible red at 655 nm wavelength with infrared laser of 830 nm.

TECH Laser is very useful for treating patients with Pain, inflammation, pressure sores and ulcers. TECH Laser has excellent wound healing properties. Useful for patients from departments of Orthopaedics, Physiotherapy, Plastic Surgery, Dermatology, Rehabilitation, and Sports Medicine. Athermic Laser Light acts on the cellular level as a bio activator. There is no heat, vibration or unpleasant electric currents.

### FEATURES

- ☀ Visible red Laser at 655 nm with output power 24 mw.
- ☀ Infrared Laser at 830 nm with output power 120 mw.
- ☀ Scanning mode: Very useful for treating larger areas such as shoulder and lumbar region without strain for the operator.
- ☀ Pulsed and continuous output.
- ☀ Meridian point detector for laser acupuncture.
- ☀ Computerised LCD screen with controls gives messages/information to the user on various adjustments.
- ☀ Automatic electronic treatment timer.
- ☀ Key operated switch prevents unauthorised operation.
- ☀ Mounted on stand and wheels for easy mobility.
- ☀ Supplied with protective goggles for the operative as well as for the patient.

### TECHNICAL SPECIFICATION

Visible Red Laser	:	Solid State 655 NM Wave Length
Power Output	:	24 MW
Infra Red	:	Solid State 830 NM Wave Length
Power Output	:	120 MW
Output	:	Cont / Pulsed

### STANDARD ACCESSORIES

- ☀ Laser protective Goggles
- ☀ Mains cable
- ☀ Meridian Point detector
- ☀ Inactive electrode