

MICROSTIM

Genius



**COMPUTERISED
MUSCLE
STIMULATOR
WITH TENS**

MICROSTIM

COMPUTERISED MUSCLE STIMULATOR WITH TENS

Genius

MICROSTIM Genius is an advanced microprocessor based electro-therapy equipment. It gives 69 programmes for MUSCLE STIMULATION.

LCD panel and a soft touch membrane keyboard enables easy setting of the parameters. The attractive compact case makes it ideal for hospitals, clinics, homes, office and travel.

SPECIAL FEATURES

MICROSTIM, Genius is a user friendly unit featuring HELP facility and therefore simple and easy to operate.

Liquid crystal display with yellowish Backlight displays the parameter chosen, the treatment time and the current in milliamps.

Compact and light weight makes the unit portable, particularly designed to be accommodated in a briefcase.

Automatic Intensity control ensures very safe operating environment.

Microstim Genius is upgradable.

TECHNICAL SPECIFICATION

Input Data

Mains Supply	:	240 V/AC/50 Hz
Power consumption	:	46 Watts
Main Fuse	:	1 amps
Weight	:	3 kgs.

FUNCTIONAL DATA

MUSCLE STIMULATOR FUNCTIONS:

- Plain galvanic
- Pulsed galvanic
- Functions (MST)
- Faradic current
- Surged Faradic
- Triangular
- Tens
- Surged Tens
- SD-Curve programs.

No. of programmes : 01-69

