

ACU TENS - 4



ADVANCED FOUR CHANNEL TENS UNIT

ACU TENS - 4


- **ACU TENS-4** is an advanced microprocessor based electro-therapy equipment. The attractive compact case makes it ideal for hospital, clinic, home, office and travel.

TECHNICAL SPECIFICATIONS

+ INPUT DATA

Mains Supply : +12v, 500mA
Power Consumption : 5 Watts

+ FUNCTIONAL DATA

Output Current : 0 to 60mA
Output Voltage : 0 to 100v (Peak to Peak)
Frequency : 1 to 200Hz
Pulse Width : 75 to 360μs
Output Waveform : Symmetric 
Treatment time : 5 - 20 Minutes
Alarm : The Buzzer alarms when the INTENSITY CONTROL is not at zero.

+ MECHANICAL DATA:

Dimension : 65 x 210 x 155 (H x D x W) in mm.
Weight : 1.26 Kg.
Operating Temperature : 0 - 70 Deg. C

CLINICAL APPLICATIONS

- ▶ Phantom Limb, Lower Extremity
- ▶ Low Back Pain
- ▶ Cervical Pain
- ▶ Cervical Osteoarthritis
- ▶ Degenerative Arthritis Cervical and Lumbar
- ▶ Lateral Rib Cage Pain
- ▶ Herpes Zoster Neuralgia
- ▶ Acute Muscle and Ligament tear - Ankle
- ▶ Post Paediatric Surgery
- ▶ Degenerative Arthritis - Knee Pain
- ▶ Recurrent Patellar Subluxation
- ▶ Low Extremity Pain
- ▶ Carpal Tunnel Syndrome
- ▶ Wrist Pain



manufactured by:

Technomed Electronics.
Plot No 41, Functional Industrial Area,
Perungudi, Chennai - 600 096
Tel : 044 - 24961854, Fax : 044 - 42020060
Email : Sales@technomedelectronics.com