




**TECHNOMED**  
ELECTRONICS



# **ELECTROSON 700**

**1 MHz Ultrasound Therapy with  
Skin Contact Detection**

**CALCIFIC BURSITIS · MYOSITIS · PLANTAR FASCIITIS · TENDINITIS**


# ELECTROSON 700

## KEY FEATURES

- Light Weight and Portable machine
- Large LCD Display
- Skin Contact Detection System
- Ergonomic Handle

## SKIN CONTACT DETECTION SYSTEM

Skin Contact Detection LED Turns Green to indicate effective skin contact and turns Red when contact is poor.



## TECHNICAL SPECIFICATIONS

Maximum Output Power	2.5 W/ sq.cm
Resolution	0.1 W/ sq.cm
Operating Frequency	1 MHz
Treatment Head Size	25 mm
Modes of Operation	Continuous and Pulsed
Pulsed Modes	1:1, 1:4, 1:7, 1:10
Treatment time	0-99 min
Operating Voltage	100-240V AC ~ 50 Hz
Power Consumption	46 Watts
Dimensions	290 x 230x 80 mm (L x D x H)
Weight	1.6 Kgs



Plot No: 41, Functional Industrial Estate,  
Perungudi, Chennai - 600 096.  
Tel: 044 - 42020060 | 044 - 43593498  
Email: sales@technomedelectronics.com

