

ACUTENS¹⁰⁰⁰

Advanced 4 Channel
TENS with NMES

Ideal for

NEUROGENIC
PAIN

NOCICEPTIVE
PAIN

MUSCLE
STRENGTHENING



TECHNOMED
ELECTRONICS

KEY FEATURES

- **4 Channels with Independent Controls**
- **Indication Based**
- **Advanced TENS machine with 10 modes**
- **Neuro Muscular Electrical Stimulation (NMES) also Included**
- **Smooth Intensity Control**
- **Easy to Use Menu Based Interface**

TECHNICAL SPECIFICATIONS

Number of Channels	4 – With Independent Controls
Output Current	0 to 99 mA
Output Voltage	0 to 80 V (Peak to Peak @ 1KΩ)
Frequency	1 to 300 Hz
Pulse Width	75 to 700 μs
Treatment time	0-99 min with Alarm
Modes	TENS – Monophasic, Biphasic Symmetric, Biphasic Asymmetric, Burst TENS, Frequency Modulation, Pulse Width Modulation, Normal, High TENS, Low TENS, Brief TENS NMES
Number of Indications	TENS - 27 , NMES -7. Total 34
Number of Modes	TENS– 10,NMES– 1. Total - 11
Operating Voltage	100-240V AC ~ 50 Hz
Power Consumption	25 Watts
Dimensions	290 x 230x 80 mm (L x D x H)
Weight	1.4 Kgs



TECHNOMED
ELECTRONICS

Plot No: 41, Functional Industrial Estate,
Perungudi, Chennai - 600 096.
Tel: 044 - 42020060 | 044 - 43593498
Email: sales@technomedelectronics.com

