

MICROSTIM

Genius



**COMPUTERISED
MUSCLE
STIMULATOR
WITH TENS**

MICROSTIM

COMPUTERISED MUSCLE STIMULATOR WITH TENS

Genius

MICROSTIM Genius is an advanced microprocessor based electro-therapy equipment. It gives 69 programmes for MUSCLE STIMULATION.

LCD panel and a soft touch membrane keyboard enables easy setting of the parameters. The attractive compact case makes it ideal for hospitals, clinics, homes, office and travel.

SPECIAL FEATURES

MICROSTIM, Genius is a user friendly unit featuring HELP facility and therefore simple and easy to operate.

✖
Liquid crystal display with yellowish Backlight displays the parameter chosen, the treatment time and the current in milliamps.

✖
Compact and light weight makes the unit portable, particularly designed to be accommodated in a briefcase.

✖
Automatic Intensity control ensures very safe operating environment.

✖
Microstim Genius is upgradable.

TECHNICAL SPECIFICATION

Input Data

Mains Supply : 240 V/AC/50 Hz

Power consumption : 46 Watts

Main Fuse : 1 amps

Weight : 3 kgs.

FUNCTIONAL DATA

MUSCLE STIMULATOR FUNCTIONS:

- ☀ Plain galvanic ☀ Pulsed galvanic ☀ Functions (MST)
- ☀ Faradic current ☀ Surged Faradic ☀ Triangular
- ☀ Tens ☀ Surged Tens ☀ SD-Curve programs.

No. of programmes : 01-69



Manufactured by:

Technomed Electronics.

Plot No41, Functional Industrial estate,
Perungudi, Chennai - 600 096.

Tel : 044-24961854, Fax: 044-42020060.

E-mail: sales@technomedelectronics.com