

TRAC

MOBILUX-5

*5 Section Motorised
Mobilization Table cum
Advanced Spinal Traction Table*



TECHNOMED
ELECTRONICS



**PRONE | SUPINE | SIDELINE | CERVICAL | PELVIC TILT | MOTORISED PELVIC WEDGE
MOTORISED HEIGHT ADJUSTMENT | ADJUSTABLE TRACTION PULLING ANGLE**

ONE TABLE FOR MOBILIZATION & TRACTION

Introducing the Mobilux-5 Trac. A versatile Mobilization table with an option to attach a powerful traction machine.

The Mobilux-5 Trac is powerful, versatile, rugged and long lasting. Built with space and cost saving in mind, the table is a great value addition for you.

SOFT LEATHERETTE

Mobilix-5 Trac comes with Soft 2-inch Foam with Leatherette for comfortable treatment.

MOTORISED ADJUSTMENTS


Motorised Height adjustments and Pelvic Wedge. Gas spring assisted Cervical and Pelvic Flexion.

LARGE BREATHER HOLE

The large, removable breather hole is designed to be ideal for all patient heights.

RUGGED DESIGN

Built to last, the Mobilux-5 Trac can support upto 150 Kgs. making it one of the most versatile.



VERSATILE FOR VARIOUS TREATMENTS

Mobilux-5 Trac is extremely resourceful to physiotherapists. Its flexible positioning adjustments help to offer a wide range of treatments.

Starting with Prone and Supine Treatments, coupled with pelvic flexion helps in offering superior tractions and treatments with less effort. The breather hole helps in positioning the patient effectively for the treatment. You can position the table for seated treatments, Knee exercises, cervical traction and much more.



RANGE OF MOTION

Mobilux's flexibility is its strength. It makes the table so versatile that its one-size fits all approach is the best addition to your clinic.

OTHER TECHNICAL SPECIFICATIONS

Overall Dimensions : 8 Ft x 2.5 Ft x 3 ft
Total Weight : 85 Kgs
Maximum Patient Weight : 150 Kgs



TECHNOMED
ELECTRONICS

Plot No: 41, Functional Industrial Estate,
Perungudi, Chennai - 600 096.

Tel: 044 - 42020060 | 044 - 43593498

Email: sales@technomedelectronics.com

